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REFERENZ

[1] Robert D. Nerenz, Haowei Song, and Ann M. Gronowski. Screening Method to Evaluate Point-of-Care Human Chorionic Gonadotropin (hCG) Devices for Susceptibility to the Hook Effect by hCG βCore Fragment: Evaluation of 11 Devices. Clinical Chemistry 60:4.667-674 (2014)

[2] Berger P, Paus E, Hemken PM, Sturgeon C, Stewart WW, Skinner JP, Harwick LC, Saldana SC, Ramsay CS, Rupprecht KR, Olsen KH, Bidart JM, Stenman UH. Candidate epitopes for measurement of hCG and related molecules: the second ISOBM TD-7workshop. Tumor Biol. (2013) 34:4033- 4057

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RÉFÉRENCES

[1] Robert D. Nerenz, Haowei Song, and Ann M. Gronowski. Screening Method to Evaluate Point-of-Care Human Chorionic Gonadotropin (hCG) Devices for Susceptibility to the Hook Effect by hCG β Core Fragment: Evaluation of 11 Devices. Clinical Chemistry 60:4.667–674 (2014)
[2] Berger P, Paus E, Hemken PM, Sturgeon C, Stewart WW, Skinner JP, Harwick LC, Saldana SC, Ramsay CS, Rupprecht KR, Olsen KH, Bidart JM, Stenman UH. Candidate epitopes for measurement of hCG and related molecules: the second ISOBM TD-7workshop. Tumor Biol. (2013) 34:4033–4057

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RIFERIMENTO

[1] Robert D. Nerenz, Haowei Song, and Ann M. Gronowski. Screening Method to Evaluate Point-of-Care Human Chorionic Gonadotropin (hCG) Devices for Susceptibility to the Hook Effect by hCG βCore Fragment: Evaluation of 11 Devices. Clinical Chemistry 60:4.667–674 (2014)
[2] Berger P, Paus E, Hemken PM, Sturgeon C, Stewart WW, Skinner JP, Harwick LC, Saldana SC, Ramsay CS, Rupprecht KR, Olsen KH, Bidart JM, Stenman UH. Candidate epitopes for measurement of hCG and related molecules: the second ISOBM TD-7workshop. Tumor Biol. (2013) 34:4033–4057

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HCG Pregnancy Test (Strip) (25 mIU/ml)

For Self-testing Use
Specimen: Urine
Format: Strip

CE 0197

INTENDED USE

HCG Pregnancy Test is a self-testing immunoassay made for the rapid, visual and qualitative determination of human chorionic gonadotropin (HCG) in urine specimen to aid in the early detection of pregnancy.

SUMMARY

Human chorionic gonadotropin (HCG) is a glycoprotein hormone secreted by the developing placenta shortly after fertilization. In normal pregnancy, HCG can be detected in serum as early as 7 days following conception. The concentration of HCG continues to rise rapidly, frequently exceeding 100 mIU/mL by the first missed menstrual period and peaking in the 100,000-200,000 mIU/mL range by 10-12 weeks into pregnancy. The appearance of HCG soon after conception and its subsequent rise in concentration during early gestational growth make it an excellent marker for the early detection of pregnancy.

PRINCIPLE

The HCG Pregnancy Test is a chromatographic immunoassay (CIA) for the rapid qualitative determination of HCG in urine. The membrane is pre-coated with anti- α HCG capture antibody on the test line region and goat anti-mouse on the control line region. During testing, the urine specimen is allowed to react with the colored conjugate (mouse anti- β HCG monoclonal antibody-colloidal gold conjugate), which has been pre-dried on the test strip. The mixture then moves upward on the membrane chromatographically by capillary action. For a positive result, a pink-colored line with the specific antibody-HCG-colored conjugate complex will form in the test line region of the membrane. Absence of this pink-colored line in the test line region indicates a negative result. Regardless of the presence of HCG, as the mixture continues to move across the membrane to the immobilized goat anti-mouse, a pink-colored line at the control line region will always appear. The presence of this pink-colored line serves as: 1) verification that sufficient volume is added. 2) that proper flow is obtained, and 3) as a control for the reagents.

COMPOSITION

The test strip contains anti- β HCG on a colloid gold particles and a combination of anti- α HCG coated on the membrane.

REAGENTS AND MATERIALS PROVIDED

1. HCG Pregnancy Test.
2. One desiccant.
3. One Instruction.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Clock or Timer.
2. Sample container.

WARNING AND PRECAUTIONS

1. Read the instruction carefully before performing this test. Pay attention to the position of the C and T line.
2. Do not use beyond the labeled expiration date.
3. The test kit should remain in the sealed pouch until use. Do not use if pouch is damaged or opened.
4. Do not reuse the test kit. Discard it in the dustbin after single use.
5. Do not touch the membrane located within the windows.
6. Do not swallow the desiccant.

STORAGE

The test strip should be stored at temperature 2-30°C, the sealed pouch for the duration of the shelf time (36months). Do the test in 1 hour when you open the pouch. DO NOT FREEZE.


SAMPLE COLLECTION

1. Collect fresh urine sample by using a disposable container which is clean and dry.
2. First morning specimens generally contain the highest concentration of HCG for early detection of pregnancy. However, any urine specimen is suitable for testing.
3. Specimens may be kept at room temperature for 8 hours. If the sample cannot be tested immediately, may store the sample at 2-8°C for 48 hours or at -20°C for a long time. Do not make the specimens repeat freeze and thaw.
4. Bring the urine sample to room temperature before testing.

TEST PROCEDURE

Allow the test strip, urine specimen and/or controls to equilibrate to room temperature (15-30°C) prior to testing. Do not open pouches until ready to perform the assay.

1. Remove the test strip from the sealed pouch and use it as soon as possible.
2. Immerse the test strip in the urine sample for at least 5 seconds with the arrow end pointing toward the urine. Do not immerse the strip above the printed MAX line.
3. Take the strip out and lay the strip flat. Or leave the test strip in the urine specimen as long as the strip is not immersed above the MAX line.
4. Wait for pink-colored lines to appear. Read result within 5 minutes. Do not read result after 5 minutes.

INTERPRETATION OF RESULTS

PREGNANT: Two distinct pink-colored lines appear, one in the test region (T) and one in the control region (C).

NOTE: The intensity of the pink color in the test region (T) may vary depending on the concentration of HCG present in the specimen. Therefore, any shade of pink color in the test region (T) should be considered positive.

NOT PREGNANT: Only one pink-colored line appears in the control region (C). No apparent pink line appears in the test region (T).

INVALID: Control Line fails to appear.

NOTE: Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test strip. If the problem persists, please contact your local distributor.

QUALITY CONTROL

A procedural control is included in the test. A colored line appearing on the control region (C) is considered an

internal positive procedural control, indicating proper performance and reactive reagents. A clear background in the membrane is considered an internal negative procedural control. If the test has been performed correctly and reagents are working properly, the background will clear to give a discernible result.

EXPECTED VALUES

Negative result will be found in the urine of healthy men and healthy non-pregnant women. However, healthy pregnant women have HCG present in their urine and serum specimens. The amount of HCG will vary greatly with gestational age and between individuals. The HCG Pregnancy Test Strip has an analytical sensitivity of 25 mIU/mL, and is capable of detecting pregnancy as early as 1 day after the first missed menses.

PERFORMANCE CHARACTERISTICS

Analytical Sensitivity: No less than 25 mIU/ml
Analytical Specificity: The test results show negative for the 500 mIU/ml hLH, 1000 mIU/ml hFSH and 1000 µIU/ml hTSH samples.
Diagnostic sensitivity and Diagnostic specificity: This HCG Pregnancy Test detects HCG at a concentration of 25 mIU/ml or greater. 900 know negative urine samples were equally divided into 6 groups. Each group of samples (150) were spiked with HCG to the concentration of 0 mIU/ml, 5 mIU/ml, 10 mIU/ml, 25 mIU/ml, 50 mIU/ml and 5 IU/ml separately, calibrated against WHO 4th international standard. Each group of sample was tested with HCG Pregnancy Test. The results from this study gave > 99% agreement with the expected results.

Result	0 mIU/ml	5 mIU/ml	10 mIU/ml	25 mIU/ml	50 mIU/ml	5 IU/ml	Total
Positive	0	0	0	150	150	150	450
Negative	150	150	150	0	0	0	450
Total	150	150	150	150	150	150	900

Diagnostic sensitivity=100% (450/450)
Diagnostic specificity=100% (450/450)
Interference Testing: The following substances were added in HCG free and HCG spiked urine samples. None of the substances at concentration tested interfered in the assay. For example:

Acetaminophen	0.2 mg/ml	Gentesic Acid	0.2 mg/ml
Acetylsalicylic Acid	0.2 mg/ml	Glucose	20 mg/ml
Ascorbic Acid	0.2 mg/ml	Hemoglobin	10 µg/ml
Atropine	0.2 mg/ml	Tetracycline	0.2 mg/ml
Caffeine	0.2 mg/ml		

LIMITATION

- Although it is not necessary to test with an early morning urine sample, excessive fluid intake should be avoided before testing. A “Not Pregnancy” result may be obtained if the urine sample is too dilute.
- The contents of this kit are for use in the qualitative detection of HCG in urine only.
- A specimen with a low level of HCG may show color development over time. If a negative result is obtained but pregnancy is suspected, another specimen should be collected after 48-72 hours and tested.
- Fertility drugs containing HCG can give misleading results (these fertility drugs are usually given by injection and testing too soon after administration may give a false “pregnancy” result).
- Other fertility therapies (such as clomiphene citrate), painkillers and hormonal contraceptives (e.g. contraceptive pill) should not affect the result.
- HCG may remain detectable for a few days to several

- weeks after delivery, spontaneous abortion, or HCG injections.
- Ectopic pregnancy, ovarian cysts, menopause, and some very rare medical conditions can give misleading results.
 - While pregnancy is the most likely reason for the presence of HCG in serum and urine, elevated HCG concentrations unrelated to pregnancy have been reported in some patients, for example, trophoblastic disease and certain nontrophoblastic neoplasms.
 - This test provides a presumptive diagnosis for pregnancy. A confirmed pregnancy diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated.
 - Alcohol may interfere the test result. It is not recommended using the test after drinking.
 - The user with eye diseases such as color blindness or color weakness may interpret the result wrongly.

REFERENCE

[1] Robert D. Nerenz, Haowei Song, and Ann M. Gronowski. Screening Method to Evaluate Point-of-Care Human Chorionic Gonadotropin (hCG) Devices for Susceptibility to the Hook Effect by hCG βCore Fragment: Evaluation of 11 Devices. Clinical Chemistry 60:4.667–674 (2014)
[2] Berger P, Paus E, Hemken PM, Sturgeon C, Stewart WW, Skinner JP, Harwick LC, Saldana SC, Ramsay CS, Rupprecht KR, Olsen KH, Bidart JM, Stenman UH. Candidate epitopes for measurement of hCG and related molecules: the second ISOBM TD-7workshop. Tumor Biol. (2013) 34:4033-4057

INDEX OF SYMBOLS

	Do not re-use		Batch code
	In vitro diagnostic medical device		Use-by date
	Store at 2-30°C		Contains sufficient for <n> tests
	Manufacturer		CE Mark
	Authorized representative in the European Union		Consult instructions for use